

# Installation Instructions for Free Standing FRP Cat Tower System

# Introduction

Thank you very much for your investment in Mason kennels. At Mason we take great pride in providing our customers with the highest quality animal enclosures combined with an enjoyable ordering experience. The following instructions will assist you with proper assembly, cleaning, and maintenance of your Mason kennels. It is important to follow these guidelines to receive the best results and maximum life from your investment.

#### **Assembly**

The following pages will show you how to assemble your new Mason kennel. Since every order is custom built to your specifications, the instructions are designed to show the various methods used to assemble our kennels and some of the demonstrations may not apply to your order. Be sure to use the enclosed scale drawings at the end of this document of the provided hardware to ensure proper identification and usage.

#### **Hardware**

Depending on your Kennel project you may have some or all the following hardware:

| Part #  | Description                               | Part # | Description                               |  |  |
|---|---|--------|---|--|--|
| 2222  | #10-32 x 5/8" TEK screw(s)                | 8081   | Triple SS clamp(s) -1" (w/stainless steel |  |  |
| 1955  | 5/16" x 1-1/2" SS carriage bolt(s)        |        | gate option)                              |  |  |
| 356   | 5/16" x 1-1/2" carriage bolt(s)           | 8082   | Double SS clamp(s) -1" (w/stainless steel |  |  |
| 1771  | 5/16" x 1- 3/4" SS carriage bolt(s)       |        | gate option)                              |  |  |
| 357   | 5/16" x 1-3/4" carriage bolt(s)           | 8083   | Single SS clamp(s) -1" (w/stainless steel |  |  |
| 424   | 5/16"-18 SS hex nut                       |        | gate option)                              |  |  |
| 423   | 5/16" -18 nut                             | 330    | Panel Clamp (w/Chain link Gate option)    |  |  |
| 5757  | #10-12 x 3/4' SS truss head screw         |        |   |  |  |
| 369   | 1/4"-20 x 1-1/2" hex head bolt            | 328    | Triple Clamp (w/Chain link Gate option)   |  |  |
|   | (w/Stainless Steel Gate option)           |        |   |  |  |
| 422   | <sup>1</sup> / <sub>4</sub> " -20 hex nut | 3118   | Single Clamp -1" (w/Stainless Steel Gate) |  |  |
|   |   |        | Option)                                   |  |  |
| 1400  | 1/4" x 1-1/2" Rawl Spike(s)               | 3119   | Double Clamp -1" (w/Stainless Steel Gate  |  |  |
|   | (anchored systems only)                   |        | option)                                   |  |  |
| 2449  | Sealant                                   | 3120   | Triple Clamp -1" (w/Stainless Steel Gate  |  |  |
|   |   |        | option)                                   |  |  |
| 1071  | Special pulleys (w/Kennel Door            |        |   |  |  |
|   | option)                                   |        | Tie Plate(s) $1-1/2$ " x length of panel  |  |  |
|   |   |        |   |  |  |
|   |   |        | "Sanislope T"                             |  |  |
| 1344  | 1/4" X 1-1/4 "Hex Head Tapcons            | 421    | 1/4" Nyloc nut                            |  |  |
| 3399  | Shelf Retaining Tab                       | 1906   | 1/4-20 x 5/8" Flat head screw             |  |  |
| Note: Wall mounting hardware is shown only for reference purposes and is not included |   |        |   |  |  |

#### Installation

## **Required installation tooling:**

- Level (4' preferred)
- 1 /2" Combination wrench
- 7/16" Combination wrench
- Hammer
- Hammer drill (anchored systems only)
- 1 /4" Masonry bit (anchored systems only)
- 5/16" drill bit (anchored systems only)
- #2 Phillips bit (2-3)

- Caulk gun
- Variable Speed Drill (cordless preferred)
- 11/64" drill bit
- Tape measure
- Chalk line
- Denatured alcohol
- Shop rags
- 5/16" Nut driver bit


#### **Optional tools**

- Vise grips 8" locking C-clamps (3 pair recommended)
- WD40 or other thread lubricant

## **Section A Panel Identification**

Refer to the floor plan supplied.


Depending on the nature of your Kennel Project you will have either letter and or number designations on the tag(s) wired to the top of the panel as shown in Figure 1.



Letter / Number designation

Figure 1

Each connection point of the panels should be aligned as shown on the floor plan. Note the difference in the alignment of the two panel connections shown in Figure 2. It will make a difference as to how your runs go together if you do not align the panels properly and in some cases they may not go together at all.



Form #159 Rev. G

# Section B Panel Assembly and Floor Layout Preparation

Sometimes the length or layout of the division panels in a kennel requires that they be made up from multiple panels. Each of the multiple panel assemblies must be connected before assembling the runs. Referring to the Floor Plan and using the Panel Identification Tag as described in Section A, make all multi-division panel connections that are in line with one another using the following steps:

- 1. Lay the two panels to be connected on the floor so that they are oriented correctly per the Floor Plan and or "Panel Identification".
- 2. Align the top surfaces and pull the mating edges together tightly. Center the "tie plate" over the seam between the two panels as shown in Figure 3. The tie plate can be applied to either side of the panels but for appearance sake you may want to apply them consistently to the same side. Hold the panels and the tie plate so that they are flush at the top edge of the division panels and tightly pulled together.
- 3. Using the tie plate as a template, drill an 11/64" hole approximately 5/8" into the panels (not all the way through) at each hole in the tie plate.
- 4. Using the #10-32 x 5/8" TEK screws provided, connect the panels by applying a screw in every hole drilled in step B-3.
- 5. Repeat steps B-1 through B-4 for all inline division panels.




Figure 3

#### **Section C**

#### **Run Assembly**

- 1. Refer to the floor plan and select a starting point for assembling the runs (preferably on the end of a bank of runs at the intersection of the back panel and outside division panel), which is best for your application.
- 2. If your system is using Silvis Seals you must apply them to the Sanislope T's now. Orient the "T" as shown in Figure 4. Apply the end blocker by peeling off the backing strip and positioning the blocker at the end of the "T".




Figure 4

3. Apply the Silvis Seal to the "T" by removing approximately 2 ft. of backing tape. Starting at the end blocker, apply the Seal to the recessed channel in the "T" being careful to orient the Seal in line with the "T". Lay the rest of the Seal in the recessed channel. Carefully pull the remaining backing tape from the "Seal" and using light pressure adhere it to the aluminum as shown in Figure 5.




Figure 5

- 4. Layout the "Sanislope T" sections approximately where they will lie on the final floor plan as shown in Figure 6. Note: be sure that the end blocker that was installed in step 2 is at the outside or gate end of the run.
- 5. **Optional**: This system is designed to function equally well whether it is fastened to the floor or not. If you do not want to fasten it the floor skip to step 7; otherwise proceed with step 6.
- 6. Using the 5/16" drill bit, pre-drill the "Sanislope T's" for the outside division panels only 3" from each end, one hole in each side, diagonal from each other. These holes will be used later to anchor the system to the floor.
- 7. Starting from the point you selected in step C-1 place the first "outside division panel" on the corresponding "Sanislope T". Place the adjacent back panel as shown in Figure 7. Note that the vertical edges of the panel are flush with each other.




Figure 7

- 8. Begin assembling this panel to the first back panel using the #10-32 x 5/8" TEK screws provided, connect the panels by applying a screw in every hole in the attaching angle to the back panel. Note: these are self-drilling, self-tapping screws and do not require a pilot hole in the division panels.
- 9. Position the next division panel in place and connect with the #10-32 x 5/8" TEK screws provided. See Figure 8.





Figure 8

- 10. Repeat steps C-1 through C-6 to assemble all the division panels to the back panels.11. Repeat steps B-1through C-6 for all other banks of runs.

# Section D Installing the floors

1. Install 2 shelf retaining tabs (#3399) to the bottom of each floor using #10-12 x .75" phillips truss head screws. USE THE HOLE NOT THE SLOT IN THE TAB. Turn the screws in until the safety tabs are pulled securely to the floor, then back the screws out approximately 1/8 of a turn. This will allow the tab to be turned while retaining a small amount of tension. See Figure 9.



2. Place the floors on the aluminum shelf support flanges in the order shown on your plans. Rotate the safety tabs to point toward the sides of the division panel and fit under the aluminum shelf support. See figure 10.




Figure 10

3. The shelves are intended to be removable for cleaning and should not be permanently fastened in place or sealed to the division panels or the back panel. To remove shelves, turn the retaining tabs one quarter turn in either direction.

# Section E Gate / Stall Front Hanging

Your system may include Chainlink or Aluminum Framed Stainless Steel gates or stall fronts or both. Use the appropriate instructions for hanging your gates and or stall fronts.

#### **Chain link Gates / Stall Fronts**

1. Attach the filler pipes to the front end of each division panel with two Square-Round clamps (it takes two halves to make up one clamp) approximately 6" from either end of the filler pipe using the 5/16" x 1-3/4" carriage bolts and 5/16" nuts provided as shown in Figure 11.




Figure 11

- 2. Attach the gate / stall front on the front of each run between the filler pipes with two clamps on each side using the 5/16" x 1-1/2" bolts and nuts provided as shown in Figure 12. Note: It takes two halves to make up one clamp. Attach gates / stall fronts to outside division panels using "panel clamps". All other gates / stall fronts will use "triple clamps". To help ensure the safety of the animals the round head of the bolt should be assembled toward the inside of the run. Note: The triple clamps attach two gates / stall fronts together. Therefore, you will have to hang two gates / stall fronts at the same time (see Figure 13).
- 3. Adjust the gate latch following the instructions in section H, "Gate / Stall Front Adjustment".





Figure 12



#### **Stainless Steel Gates / Stall Fronts**

- 1. Attach the gate/stall front on the front of each run between the division panels with three clamps on each side using the 5/16" x 1-3/4" bolts and nuts provided as shown in Figure 14. Note: It takes two halves to make up one clamp. Attach gates/stall fronts to outside division panels using "square 1" panel clamps". All other gates/stall fronts will use "square 1" triple clamps". To help ensure the safety of the animals the round head of the bolt should be assembled toward the inside of the run. Note: The triple clamps attach two gates/stall fronts together. Therefore, you will have to hang two gates/stall fronts at the same time (see Figure 15).
- 2. Adjust the gate latch following the instructions in section H, "Gate / Stall Front Adjustment".

|   | 2 |   |     |
|---|---|---|-----|
| 4 | 5 |   |     |
| 1 |   |   |     |
|   |   |   | 2 2 |
| 0 |   | 1 | 5 5 |

| Parts List       |      |                      |  |  |  |
|------------------|------|----------------------|--|--|--|
| ITEM PART NUMBER |      | DESCRIPTION          |  |  |  |
| 1                | 424  | 5/16-18 SS nut       |  |  |  |
| 2                | 1771 | 5/16-18 x 1.75" SS   |  |  |  |
|                  |      | carriage bolt        |  |  |  |
| 3                | 8081 | 1in to 0.75in_Triple |  |  |  |
|                  |      | Gate Clamp           |  |  |  |
| 4                | 8082 | 1in to               |  |  |  |
|                  |      | 0.75in_Double Gate   |  |  |  |
|                  |      | Clamp                |  |  |  |
| 5                | 8083 | 1in to 0.75in_Single |  |  |  |
|                  |      | Gate Clamp           |  |  |  |

Figure 14




Figure 15

This Mason system is designed to work equally well whether it is anchored to the floor or not. If you desire to anchor it to your floor proceed with Section F otherwise skip to Section G.


#### **Section F**

#### **Anchor System to Floor**

**Important:** Before anchoring your system to the floor verify the entire installation is where you want it, square to the lines marked in step C-1 and leveled. Once you have verified this, you can anchor the "Sanislope T's" on the outside division panels to the floor using the following steps:

- 1. Using the holes you drilled in step C-3 as a template, drill a 1/4" hole using a hammer drill, with a 1/4" masonry bit to a minimum of 1-1/4" deep.
- 2. Insert a 1/4" x 1-1/2" Rawl Spike provided into the hole and drive it in with a hammer until the "Sanislope T" is tight against the floor and the head of the Rawl Spike is tight against the "Sanislope T".
- 3. Repeat steps F-1 and F-3 for all hole locations shown in Figure 16.
- 4. Repeat steps F-1 through F-4 for all outside division panels.

#### DETAILS OF ALUMINUM FRAMED SANI-KENNEL PANEL



# Section G Sealing

- 6. Clean along the edge of all of the Division and Back panels where they come in contact with the concrete using denatured alcohol applied to shop rag.
- 7. Apply a thin bead of sealant provided along the seam created between the panels and the concrete as shown in Figure 16.

#### **Section H**

## **Gate / Stall Front Adjustment**

During shipment and installation, gate assemblies may be forced out of adjustment. If your gate does not latch automatically when pushed closed, the latch can be adjusted by following these instructions.

#### <u>Chain link Gate / Stall Front</u> (Reference Figure 17)

During the following adjustment process you may need to move the wire mesh slightly. If necessary, use a small block of wood and hammer to tap the wire out of the way.

- 1. Loosen the latch catch bolt just enough to allow the latch catch to slide up and down.
- 2. Tap the latch catch up or down until the inside latch bar handle rests on the horizontal brace pipe and the latch bar rests in the bottom of the latch catch.
- 3. Check to insure the pendant swings freely. If it catches on the latch catch bolt or you've run out of adjustment in the latch catch slot you may have to raise the gate inside of its frame.
- 4. If the gate needs to be moved, loosen the top and bottom hinge bolts and raise (or lower) both the gate and the top hinge as needed then tighten both hinge bolts.
- 5. Re-adjust the latch catch.
- 6. Tighten latch catch bolt.




Figure 17

# **Stainless Steel Gate / Stall Front** (Reference Figure 18)

- 1. Loosen the latch catch bolts just enough to allow the latch catch to slide up and down.
- 2. Tap the latch catch up or down until the inside latch bar handle rests on the inside of the keeper and the latch bar rests in the bottom of the latch catch.
- 3. Tighten latch catch bolts.






Figure 18

#### **Maintenance**

Over time, your Mason kennels might require adjustments, lubrication, or replacement parts in order to remain in top working condition. We recommend a yearly maintenance schedule to lubricate door hinges, adjust gate locks and any other moving parts. A standard, industrial grease works well. If your kennels should need any replacement parts, our professional sales engineers will be happy to review your original order and assist you.





Form #159 Rev. G