

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

DANGER. Corrosive. Causes irreversible eye damage. Harmful if swallowed. Do not get in eyes or on clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove contaminated clothing before reuse.

FIRST AID

If in Eyes: Hold eye open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses, if present, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

If Swallowed: Call poison control center or doctor for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by the poison control center or doctor. Do not give anything by mouth to an unconscious person.

Have the product container or label with you when calling poison control center or doctor, or going for treatment.

NOTE TO PHYSICIAN

Probable mucosal damage may contraindicate the use of gastric lavage.

ENVIRONMENTAL HAZARDS

This product is toxic to fish and other aquatic organisms.

PHYSICAL/CHEMICAL HAZARDS

Keep away from combustible, organic or other readily oxidizable materials. Do not mix with any other chemicals.

VistaTab™ Dental Waterline Cleaner Tablets

by Hu-Friedy

- Cleans and Protects Dental Unit Waterlines
- For the cleaning and control of microbial contamination in dental unit waterlines

Active Ingredients:

Sodium Chlorite 4.0%
 Sodium Dichloroisocyanurate dihydrate. 1.0%
 Other Ingredients 95.0%
 TOTAL..... 100.0%

EPA Reg. No. 70060-25-75372
EPA Est. No. 74245-AZ-1

*Sodium chlorite in this product is from Spain.

KEEP OUT OF REACH OF CHILDREN
DANGER

This box contains: 15 - 1.7 gram (0.06 ounces) tablets
Net Weight 25.5 grams (0.90 ounces)

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

INITIAL TREATMENT:

1. Add two (2) 1.7 gram VistaTab™ Dental Waterline Cleaner Tablets for every 750 mL of potable water. Allow the tablets to dissolve.
2. Flow about half of the treated water through the dental water lines and let the treated water sit within the lines for 5 minutes. Purge the remaining water through the lines.
3. Repeat steps 1 and 2.
4. Flush at least 500 mL of potable water through the lines to remove any residual chemical.

ROUTINE TREATMENT:

1. Every week or as needed*, add one (1) 1.7 gram VistaTab™ Dental Waterline Cleaner T tablet for every 750 mL of potable water. Allow the tablet to dissolve.
2. Flow about half of the treated water through the dental waterlines and let the treated water sit within the lines for 5 minutes. Purge the remaining water through the lines.
3. Flush at least 500 mL of potable water through the lines to remove any residual chemical.

Note: Periodic cleaning may be done on any day of the week.

*Treatment should be done if bacteria count is greater than 500 CFU/mL as recommended by the CDC.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

PESTICIDE STORAGE: Store this product in a cool, dry area away from direct sunlight and heat. Keep product from contacting air prior to use to avoid deterioration.

PESTICIDE DISPOSAL: Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture, or rinsate is a violation of Federal law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance.

CONTAINER HANDLING AND DISPOSAL: Non-refillable. Do not reuse or refill this packet. Offer for recycling, if available. If recycling is not available, dispose empty packet in sanitary landfill or by incineration. If burned stay out of smoke.



Manufactured for:
Hu-Friedy Mfg. Co., LLC
3232 N. Rockwell St. | Chicago, IL 60618 | USA
1-800-Hu-Friedy or 1-773-975-6100

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
©2022 Hu-Friedy Mfg. Co., LLC. All rights reserved.

IMS-1451

