

Midmark®

Instrument processing

Education Options

Thank you for your interest in learning more about Infection Prevention. At Midmark, we understand the importance of proper Clinical Education and the Instrument Processing procedures that need to take place to ensure better care for your patients. **Below you will find an overview of our education options available to you and your team!**

COMPLIMENTARY OPTIONS

- Education introduction letter with support links and reference guides.
- Scheduled discovery call designed to review your current processes so we may tailor our instrument processing education services to your team(s).

VIRTUAL AND ON-SITE OPTIONS

To be quoted after your Discovery Call

1-hour Virtual Instrument Processing Basics

- Choosing the correct cycle for your instruments properly using manufacturers' instructions for use.
- Instrument processing workflow from point-of-use to sterile storage.
- Monitoring your sterilizer and instruments based on current best practices.
- Performing maintenance on your sterilizer to ensure consistent performance.
- Documentation requirements and recommendations for a compliant sterile processing area.

4-HOUR ON-SITE INSTRUMENT PROCESSING BASICS

- Engage in practical demonstrations to apply virtual content hands-on.
- Implement learned principles through a 4-hour session.

ADDITIONAL RESOURCES INCLUDED WITH VIRTUAL OR ON-SITE OPTIONS

- eBook with support documents and competency checklist to assist with standardization.
- Certificates of completion for all attendees.
- Two-week follow-up call to seek additional guidance and share your experiences.