



# Midmark HEINE® Ophthalmoscopes

## COMPARISON CHART



PRODUCT FEATURES	HEINE BETA® 200	HEINE K180™
		
<b>Performance</b>		
Light Source	LED	LED
Color Temperature	LED: 3,300 K	LED: 3,300 K
CRI	Typ. 95 R9 typ. 93	Typ. 95 R9 typ. 93
Brightness Control	Continuous between 100 - 3% with one-finger operation	Continuous between 100 - 3% with one-finger operation
Optics	Aspherical Optical System (AOS) elliptical beam, reflex free	Precision spherical optics with proprietary glare-mitigating coating
Magnification	15x	15x
Field of View	5-10°	5-10°
<b>Durability</b>		
Sealing	Dust-resistant	Dust-resistant
Frame	Aluminum	Polycarbonate
Window	Scratch-resistant glass	Scratch-resistant glass
Light Source Life Expectancy	LED: typ. 100,000 hrs	LED: typ. 100,000 hrs
<b>Function</b>		
Lenses/Rekoss Disc	27	27
Range of Lenses	27 single lenses +1D steps from 1-10, 15, 20, 40 -1D steps 1-10, 15, 20, 25, 35	+1D steps from 1-10, 15, 20, 40 -1D steps 1-10, 15, 20, 25, 35
Apertures	5 - large spot, small spot, hemislot, slit, fixation star with grid	3 - large spot, small spot, slit
Filters	Red-free with all apertures, cobalt blue	Red-free, cobalt blue

For product inquiries, please contact the Midmark Medical Customer Experience Department at 844-856-1233, option 1.  
For more information or a demonstration, contact: Midmark Corporation, Versailles, OH.  
Call: 1-800-MIDMARK or visit: [midmark.com](http://midmark.com)  
HEINE, BETA and K180 are trademarks of HEINE Optotechnik GmbH + Co. KG.

MKT-00155

