

# Explore the next generation of Midmark® sterilizers



## Designed for better:

- DURABILITY
- WORKFLOW
- COMPLIANCE



## Increase Longevity

**25,000 CYCLES**  
Gain 2.5-times longer device life than previous models.

## Reduce Lifetime Cost

**UP TO \$11,800 SAVINGS**  
Save thousands compared to competitive models with a sterilizer designed to require less frequent routine technician-performed service.

## Enhance Workflow



**3 LANGUAGES**  
Select English, French or Spanish for your preferred display language.\*

\*French and Spanish available with a future software update



**LIGHTER, LARGER TRAYS**  
Increase usable space with a larger but lighter tray than previous models.




**5-INCH TOUCHSCREEN**  
Navigate setup with a clear, graphical display you can utilize while wearing gloves.


## Improve Efficiency

### AUDIO + VISUAL INDICATORS


Determine cycle status—cycle in progress, cycle complete and cycle error—from a distance with a progressing color-coded LED light bar, large cycle-countdown clock and audible cues.



**CYCLE IN PROGRESS**



**CYCLE COMPLETE**



**CYCLE ERROR**



**TIME-SAVING DOCUMENTATION**  
Help save time, simplify teammate training and build instrument processing confidence with automated mechanical monitoring and chemical indicator reporting.



**OPTIONAL AUTOFILL + AUTO-DRAIN**  
Get fresh water for every cycle without burdening your team.

## Reduce Stress



**ONSCREEN GUIDANCE**  
Improve safety with touchscreen prompts that lead your team through a standardized process.



**ON-DEVICE TUTORIALS**  
Reduce the burden of training with step-by-step tutorials.



**STREAMLINED COMPLIANCE RECORDKEEPING**  
Stay audit-ready with device reminders and notifications, user authentication, unlimited storage of routine care events for the life of the sterilizer, and automated cycle recordkeeping.



### MORE CYCLE OPTIONS

Simplify clinical workflow with dynamic-air-removal steam sterilizers that provide unique preset cycles to meet both the CDC and AAMI sterilization cycle recommendations for wrapped instruments.



**WRAPPED**  
270°F, 4 MIN



**DELICATE**  
250°F, 30 MIN

**NEW**



**WRAPPED**  
275°F, 3 MIN



**UNWRAPPED**  
270°F, 3 MIN



**5 CUSTOM CYCLES\*\***  
250°F–275°F, 3 MIN–45 MIN

\*\*These cycles are not FDA-cleared, and validation of the sterility of items processed using a custom cycle is the responsibility of the user.

**Simplify sterilization and get your team back to the work they love. Choose technology designed with them in mind.**

[LEARN MORE >](#)



© 2024 Midmark Corporation, Versailles, Ohio USA  
MKF00573

