



### SOPRO® INTRAORAL CAMERAS BY ACTEON

(Compatible with Midmark Procenter and Elevance® Delivery Units)

Easily capture images without blurring when you use a SOPRO intraoral camera with the SOPRO Touch sensor. SOPRO cameras are lightweight, easy to clean and compatible with both analog and digital (USB 2.0) uses. **Note:** SOPRO cameras require a docking station. Using a deluxe rheostat is recommended.



#### 617 CAMERA

- Focus light with minimal shadowing using the latest generation of LED lights with condensers.
- Improve exploration of the distal areas with a wide 105-degree view.
- Improve patient comfort with a curved design that follows the natural shape of the mouth.
- Fix your portrait down to a single tooth with a fixed focus and large depth of field.
- Avoid distortion and produce a high-quality image with an aspheric lens.



#### 717 FIRST CAMERA

- See details invisible to the naked eye with magnification up to 115 times in Macro mode without any additional equipment.
- Achieve perfect illumination with 8 LED lights.
- Capture images in 4 magnification modes: Macro (115x magnification), 0–5 mm, intraoral 3–30 mm and extraoral 30 mm+.
- Gain better access with a sleek design and extremely thin distal part.
- Achieve exceptional image quality with the digital signal processor.



#### SOPROLIFE® CAMERA

- Identify decay and provide minimally invasive treatment with the clarity of 3 visual modes (Daylight, Diagnostic and Treatment).
- Communicate more effectively with your patients using the white light images of Daylight Mode.
- Identify the development of occlusal and/or proximal carious lesions earlier using Diagnostic Mode.
- Visually differentiate between infected and healthy tissue in the excavated site using Treatment Mode.
- See details invisible to the naked eye with magnification up to 115 times in Macro mode.



#### SOPROCARE® CAMERA

- Increase case acceptance by easily highlighting caries, plaque, calculus and gingival inflammation with this fluorescence-based camera.
- Capture the patient's portrait, single tooth, quadrant, arch or smile images using Daylight Mode.
- Detect caries clearly with red coloration using Caries Mode.
- Highlight gingival inflammation, plaque and calculus using Perio Mode.

## BIEN-AIR® ELECTRIC HANDPIECES

(Compatible with Midmark Procenter and Elevance® Delivery Units)



### MX2 PLUS MICROMOTOR

- 100-40,000 RPM
- Torque: 3.5 ncm
- Weight: 79 grams (0.17 lb)
- Max sound level: 51 dBA
- Autoclavable
- Lifetime lubricated ball bearings
- Medical grade stainless steel
- 3-year warranty



### NOVA® 1:5 HIGH-SPEED CONTRA ANGLE

- 500-200,000 RPM
- Friction grip chuck
- Small tapered head
- 6.3 mm head diameter
- 13.7 mm height
- Stainless steel body
- Shockproof head
- Push-button with CoolTouch+™ heat-arresting technology
- Autoclavable anti-retraction valve
- E-type attachment (compatible with most motor brands)
- 2-year warranty



### NOVA® 1:1 LOW-SPEED CONTRA ANGLE

- 100-40,000 RPM
- Latch chuck
- Small tapered head
- 6.3 mm head diameter
- 12.9 mm height
- Stainless steel body
- Shockproof head
- Push-button with CoolTouch+™ heat-arresting technology
- Autoclavable anti-retraction valve
- E-type attachment (compatible with most motor brands)
- 2-year warranty



### STRAIGHT 1:1 NOSECONE

- Straight PM 1:1
- External water spray
- E-type attachment (compatible with most motor brands)
- 2.35 mm bur shank diameter
- Autoclavable
- 2-year warranty

## SCALERS

(Compatible with Midmark Procenter, Elevance, and Standard L/R and Console Asepsis 21® Delivery Units)



### DENTSPLY® CAVITRON®

- Enable free flowing movement of an ultrasonic insert within the oral cavity with a fully rotating handpiece.
- Increase quality and performance with Cavitron Magnetostrictive Technology.
- Maintain tip stroke with automatic system power adjustment, compensating for load and conditions.
- Receive 30 KHz insert included.
- Control on/off with rheostat.



### ACTEON® PURE NEWTRON® SLIM SCALER

- Ensure the proper tip and power setting with the color-coding system.
- Regulate frequency and power automatically with the "cruise control" setting.
- Control on/off through rheostat.
- Receive three included tips (#1, 10Z, H3).



### ACTEON MINI LED CURING LIGHT

- Combine power, efficiency and speed with an LED curing light.
- Utilize a wide light spectrum (420-480 nm).
- Choose between three modes of cure (fast, pulse and ramp).
- Receive included light shield and 7.5 mm autoclavable light guide.
- Get a stylish handpiece that is compact, lightweight and well balanced.
- Simplify safety with shock resistance, silent operation and easy cleaning.

## CURING LIGHT

(Compatible with Midmark Procenter and Elevance Delivery Units)

Acteon, Bien-Air, NOVA, CoolTouch+, Newtron, Dentsply, Cavitron, SOPRO, SoprolIFE and SoprocARE are not trademarks of Midmark Corporation and are the sole property of their respective owners.

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: [midmark.com/ISO](http://midmark.com/ISO)

For more information, contact your Midmark dealer or call: 1.800.MIDMARK Outside the USA call: 1.937.526.3662 or visit our website: [midmark.com](http://midmark.com)

© 2022 Midmark Corporation, Miamisburg, Ohio USA. Products subject to improvement changes without notice. Litho in USA.

007-10415-00 Rev B1 (3/22)