



  
midmark®

MEDICAL

# 640 Pediatric Examination Table



# Safe, secure care for **your** **most precious** **patients**



When caring for infants, moving the patient from one station to the next can create unnecessary discomfort and risk. That's why we designed the Midmark® 640 Pediatric Examination Table to weigh, measure and examine infants in one safe, secure location. Because at Midmark, we believe that who's using our equipment is every bit as important as how and why it's being used.



# Weigh, measure and examine infants **safely** and **securely**





















The Midmark® 640 is designed with features that can help you weigh, measure and examine infant patients safely, accurately and efficiently. Press the Zero/On button and within seconds the scale display will read, "Scale ready." Place the infant on the scale and the weight will automatically display. When you want to measure the patient, simply pull out the infantometer, then slide it back into place when finished.



- 01 **Reweigh** allows you to recapture weight with the press of a button—no need to pick the infant back up.
- 02 A **surface height of 36 inches** provides easy access and reduces the need for repeated bending during the exam.
- 03 **Upholstered end panels, beveled sides** and **flat surfaces with no sharp edges** help ensure patient safety.
- 04 For storage needs, there is an **oversized drawer** as well as a **bulk storage area** designed for keeping diapers and equipment.
- 05 A **hidden paper roll holder** with a paper tear strap at one end adds convenience.
- 06 The **built-in scale** measures in pounds or kilograms with unit lock to help ensure consistency and automatically shuts off after 30 seconds to conserve energy.

# COLOR OPTIONS

## STANDARD UPHOLSTERY COLORS

	Mist 813
	Stone 814
	Restful Path 858
	Iron Ore 851
	Lunar Gray 845
	Obsidian 857
	Oasis 843
	Soothing Blue 855
	Dream 847
	Mineral 844
	Healing Waters 848
	Shaded Garden 853
	Citrus 849
	Sandy Retreat 850
	Curative Copper 846
	Cranberry 859
	Dark Linen 856
	Latte 860
	Robust Brown 852
	Deep Earth 854

## ULTRAFREE UPHOLSTERY COLORS

	UltraFree Stone 865
	UltraFree Restful Path 868
	UltraFree Obsidian 867
	UltraFree Harbor 861
	UltraFree Wheatgrass 862
	UltraFree Spice 863
	UltraFree Cranberry 869
	UltraFree Dark Linen 866
	UltraFree Latte 870
	UltraFree Branch 864

# SPECIFICATIONS

## 640

Patient Weight Capacity:  
4.5 lb (2 kg) - 40 lb (18 kg)

Scale Resoluton:  
0.5 oz for pounds-ounces  
0.02 lb for pounds  
0.01 kg for kilograms

Scale Accuracy:  
< 22 lb (10 kg) +/- 1.0 oz (28 g)  
≥ 22 lb (10 kg) +/- 2.0 oz (56 g)

Built-In Measurement Arm offers  
a range of measurement 15-36"  
(38.1-91.4 cm)

Working Height:  
35.8" (90.8 cm)

Table Top Length:  
41.9" (106.4 cm)

Table Top Depth:  
20.8" (52.7 cm)

Paper Roll Holder:  
Upper portion of table  
accommodates one 14"  
(35.6 cm) long, 3" (7.6 cm)  
diameter roll

Drawer:  
33" x 15" x 4.125"  
(83.8 cm x 38.1 cm x 10.5 cm)

Drawer is constructed of high  
impact polystyrene with rounded  
corners for easy cleaning and is  
supported by side-mounted ball  
bearing steel slides.

Bulk Storage:  
33.5" x 17" x 18.38"  
(85 cm x 43.2 cm x 46.7 cm)

## SAFETY STANDARDS

640 table meets 60601  
medical safety standard for  
US and Canada.

## OPTIONS

Door/Drawer Locks



Designing better care.®



CARB 93129.2 Phase 2 Compliant and TSCA Title VI Compliant

The color examples shown are the best representation of the original material. Actual color may vary slightly. We strongly recommend that you contact Midmark Customer Experience at 1.800.MIDMARK to request a sample before placing your order.

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: [midmark.com/ISO](http://midmark.com/ISO)

For more information, contact your Midmark dealer or call: 1.800.MIDMARK

Outside the USA call: 1.937.526.3662 or visit our website: [midmark.com](http://midmark.com)

© 2021 Midmark Corporation, Miamisburg, Ohio USA  
Products subject to improvement changes without notice.  
Litho in USA.

007-10405-00 Rev. A1 (11/21)

