

There's an easier way to transfer patient data.

Streamline workflow and save time in preserving your anesthetic records.



OPTIMIZE YOUR MIDMARK MONITOR

Improve data-capture efficiency and accuracy by connecting your Midmark Multiparameter Monitor with your practice management software. The Midmark Anesthetic Record Interface streamlines the digital transfer of anesthetic records from the monitor to the patient's electronic medical record, enabling time savings on manual data transfer.

Transfer files following each procedure or as a batch at the end of the day.

The interface currently works with the following monitors and practice management software:

Midmark Multiparameter Monitors:


8019-021 through 8019-032
8020-001 through 8020-008

Practice Management Software*:

- Avimark®
- Cornerstone® (v6.0+)
- ImproMed® Advantage
- ImproMed® Advantage+
- ImproMed® (v4.6+)
- ezyVet®
- Pulse™ (formerly eVetPractice)


Midmark is continuously building interfaces with additional Practice Management Software! If you would like to subscribe to the Midmark Anesthetic Record Interface, but your clinic's Practice Management Software is not listed, please contact Midmark at midmarksubscription@midmark.com

TRANSFER PATIENT DATA INSTANTLY!



Automatically Captured Data:

- Raw vital signs data
- Anesthetic record summary
- ECG waveforms
- Individual parameter graphs



SUBSCRIBE TODAY!

Enjoy special introductory savings on your first-year subscription.

Visit midmark.com/ARI or scan the QR code to secure your annual software subscription and schedule installation.



License fees include the Midmark Anesthetic Record Interface, remote software installation and ongoing user support. One license is required per facility, regardless of the number of monitors used. There are no additional per-user fees.